

MAGNETOVET 200

MAGNETOTHERAPY

TECHNICAL FEATURES

Display

Size 2.6", backlit

Channels

2

Intensity

200 Gauss average for channel with flex solenoids
400 total average Gauss (for the 2 channels)

Frequency

From 5 to 200Hz

Duty Cycle

From 5% to 50% in steps of 5

Power supply

- Mains
- Battery pack (optional)

Size


100x160x35 mm

Weight

440 gr








EQUIPMENT

- 1 Bag
- 1 MAGNETOVET 200 unit
- 2 Flexible solenoids 300x100 mm
- 1 Power supply unit
- 1 Operating manual
- 1 User's guide



CE



-  LAST 10
-  MULTI USER
-  FAVORITES
-  PROGRAMMABLE
-  STIM LOCK
-  BATTERY PACK (optional)
-  POWER SUPPLY

MAGNETOVET 200

MAGNETOTHERAPY

Rev. En_09_2019

PROGRAMS LIST

DOG / CAT

Fractures

Pseudoarthrosis

Wounds / Ulcers

Hip osteoarthritis

Elbow osteoarthritis

Knee osteoarthritis

Shoulder osteoarthritis

Spondylarthrosis

Luxation

Cicatrization

Edema

Peripheral paralysis

Inflammation pain

HORSE

Osteoarthritis

Wounds / Ulcers

Discopathy

Edema