

# MAGNETOVET 4000

## MAGNETOTHERAPY

### TECHNICAL FEATURES

#### Display

Size 2.6", backlit

#### Channels

4

#### Intensity

200 Gauss average for channel with flex solenoids  
800 total average Gauss (for the 4 channels)

#### Frequency

From 5 to 200Hz

#### Duty Cycle

From 5% to 50% in steps of 5

#### Power supply

- Mains
- Battery pack (optional)

#### Size


170x220x60 mm

#### Weight

1370 gr

### EQUIPMENT

- 1 Bag
- 1 MAGNETOVET 4000 unit
- 4 Flexible solenoids 300x100 mm
- 1 Power supply unit
- 1 Operating manual
- 1 User's guide



- LAST 10
- MULTI USER
- FAVORITES
- PROGRAMMABLE
- STIM LOCK
- 2+2 FUNCTION
- BATTERY PACK (optional)
- POWER SUPPLY

# MAGNETOVET 4000

## MAGNETOTHERAPY

Rev. En\_09\_2019

### PROGRAMS LIST

#### DOG / CAT

Fractures

Pseudoarthrosis

Osteoarthrosis

Wounds / Ulcers

Hip osteoarthrosis

Elbow osteoarthrosis

Knee osteoarthrosis

Shoulder osteoarthrosis

Spondylarthrosis

Discopathy

Luxation

Arthritis

Arthrosis

Delayed fracture consolidation

Cicatrization

Edema

Muscle relaxation

Urinary retention

Smooth muscles

Peripheral paralysis

Inflammation pain

#### HORSE

Fractures

Wounds / Ulcers

Elbow osteoarthrosis

Shoulder osteoarthrosis

Arthritis

Arthrosis

Delayed fracture consolidation