

STIMVET 200

ELECTROTHERAPY 4 CHANNELS

TECHNICAL FEATURES

Display

Size 2.6", backlit

Channels

4 independent (8 electrodes)

Frequency

1-150Hz

Pulse amplitude

50-450µs

Intensity

80mA per channel

Power supply

Battery pack

Size


100x160x35 mm

Weight


440 gr

EQUIPMENT

- 1 Bag
- 1 STIMVET 200 unit
- 2 Grey cables for microcurrents
- 4 Colored cables for electrode connection
- 4 Silicon square electrodes
- 4 Silicon rectangular electrodes
- 4 Elastic bands
- 1 260 ml gel bottle
- 1 Power supply unit
- 1 Operating manual
- 1 User's guide



CE



STIMVET 200

ELECTROTHERAPY 4 CHANNELS

Rev. En_09_2019

PROGRAMS LIST

DOG

EMS

Recovery after ACL surgery

Contractures

Forelimb muscle atrophy prevention

Hindlimb muscle atrophy prevention

Trunk muscle atrophy prevention

Forelimb muscle reinforcement

Hindlimb muscle reinforcement

Trunk muscle reinforcement

Wobbler syndrome

TENS

Endorphinic TENS

Conventional antalgic TENS

Anti-accommodation modulated TENS

Post-surgical pain

Fracture pain

Osteoarthritis

Wobbler syndrome

Discopathy

Spondylarthritis

Lumbosacral stenosis

SPORT

Forelimbs maximum strength

Hindlimbs maximum strength

Forelimbs endurance strength

Hindlimbs endurance strength

Forelimbs explosive strength

Hindlimbs explosive strength

Decontracting

Cooling down

MICROCURRENTS

Osteogenesis

Wound cicatrization

Cartilage

CAT

EMS

Recovery after ACL surgery

Contractures

Forelimb muscle atrophy prevention

Hindlimb muscle atrophy prevention

Trunk muscle atrophy prevention

Forelimb muscle reinforcement

Hindlimb muscle reinforcement

Trunk muscle reinforcement

TENS

Endorphinic TENS

Conventional antalgic TENS

Anti-accommodation modulated TENS

Post-surgical pain

Fracture pain

Osteoarthritis

Discopathy

Spondylarthritis

Lumbosacral stenosis

MICROCURRENTS

Osteogenesis

Wound cicatrization

Cartilage