

# ULTRASOUNDVET 200

## ULTRASOUND THERAPY

### TECHNICAL FEATURES

#### Display

Size 2.6", backlit

#### Uscite

1 with Ø mm 42 mm handpiece

#### Frequency

1MHz

#### Power

1,4W/cm<sup>2</sup>, ±20%

#### Duty Cycle

From 10% to 100% in step of 10%

#### Emission

Continuous and pulsed

#### Mode

Direct contact and immersion

#### Power supply

- Mains
- Battery pack (optional)

#### Size


mm 100x160x35

#### Weight

440 gr




### EQUIPMENT

- 1 Bag
- 1 ULTRASOUNDVET 200 unit
- 1 Ø mm 42 handpiece - IP 68
- 1 260 ml gel bottle
- 1 Power supply unit
- 1 Operating manual
- 1 User's guide



CE



-  PROGRAMMABLE
-  BATTERY PACK (optional)
-  MAIN SUPPLY

# ULTRASOUNDVET 200

## ULTRASOUND THERAPY

Rev. En\_09\_2019

### PROGRAMS LIST

#### DOG (weighing more than 15 kg)

Hindlimb muscle contracture

---

Chronic osteoarthritis

---

Scar adhesions

#### HORSE

Deep tendinitis phase 1

---

Deep tendinitis phase 2

---

Muscle contractures