

ULTRASOUNDVET 4000

ULTRASOUND THERAPY

TECHNICAL FEATURES

Display

Size 2.6", backlit

Outputs

1 with Ø 42 handpiece

Frequency

1MHz and 3 MHz

Power

3W/cm², ±20%

Duty Cycle

From 10% to 100% in step of 10%

Emission

Continuous and pulsed

Mode

Direct contact and immersion

Power supply

- Mains
- Battery pack (optional)

Size


mm 170x220x60




Weight

1370 gr

EQUIPMENT

- 1 Bag
- 1 ULTRASOUNDVET 4000 unit
- 1 Ø mm 42 handpiece - IP 68
- 1 260 ml gel bottle
- 1 Power supply unit
- 1 Operating manual
- 1 User's guide



-  PROGRAMMABLE
-  BATTERY PACK (optional)
-  MAIN SUPPLY

ULTRASOUNDVET 4000

ULTRASOUND THERAPY

Rev. En_09_2019

PROGRAMS LIST

DOG

Chronic osteoarthritis

Reacutized osteoarthritis

Superficial tendinitis phase 1

Superficial tendinitis phase 2

Deep tendinitis phase 1

Deep tendinitis phase 2

Bursitis

Scar adhesions

Muscle contractures

CAT

Chronic osteoarthritis

Reacutized osteoarthritis

Scar adhesions

Muscle contractures

HORSE

Chronic osteoarthritis

Reacutized osteoarthritis

Superficial tendinitis phase 1

Superficial tendinitis phase 2

Deep tendinitis phase 1

Deep tendinitis phase 2

Bursitis

Scar adhesions

Muscle contractures