

# LASERVET 8.0

## LASER THERAPY

LASERVET 8.0 IS A PORTABLE AND HANDHELD DEVICE FOR HIGH PULSED POWER LASER THERAPY THAT ALLOWS TO PERFORM VETERINARY TREATMENTS. THE ACTION OF LASER IS AIMED AT TISSUE REGENERATION, POST-TRAUMATIC AND ANTI-INFLAMMATORY TREATMENT AND REDUCTION OF BOTH CHRONIC AND ACUTE PAIN.


SPECIFIC PROGRAMS FOR DOGS, CATS, HORSES, SMALL MAMMALS, REPTILES AND BIRDS

### TECHNICAL FEATURES

<b>Display</b>
Backlit, 2.6"
<b>Laser source</b>
Class 4 GaAlAs diode
<b>Peak power</b>
8W $\pm$ 20%
<b>Preset for wave length</b>
808nm 980nm and 1064nm (optional)
<b>Issue</b>
Pulsed and continuous
<b>Target area dimensions</b>
11 mm optics 30 mm and 60 mm optics (optional)
<b>Frequency</b>
Up to 10000Hz
<b>Power supply</b>
Mains and rechargeable battery
<b>Size</b>
100x160x35 mm
<b>Weight</b>
440 gr

### STANDARD EQUIPMENT

- 1 Carrying case
- 1 LASERVET 8.0 device
- 1 Laser handpiece 808nm with 11 mm optics
- 2 Protection glasses
- 1 Power supply unit
- 1 Operating manual









### PURCHASE OPTIONS

- Cod. G5994 Handpiece MP 980nm power 8W
- Cod. G5996 Handpiece 1064nm power 3W
- Cod. G6089 Optics 60 mm
- Cod. G6092 Optics 30 mm

CE



-  LAST 10
-  MULTI USER
-  PROGRAMMABLE
-  RENTAL OPTION
-  BATTERY PACK
-  MAINS POWER SUPPLY

# LASERVET 8.0

## LASER THERAPY

**MORE THAN 130 PROGRAMS FOR DOGS, CATS, HORSES, SMALL MAMMALS, REPTILES AND BIRDS. TREATMENTS INCLUDE PROTOCOLS FOR:**

**Muscle-skeletal injuries** such as arthrosis, tendinopathies, contusions, bursitis, sprains, muscle strains, fractures, dysplasia, neuropathies etc.

**Pathologies of the integumentary system** such as skin wounds, granulomas, sacculitis, perianal fistulas, pododermatitis, abrasions, etc.

**Pathologies of other systems** such as otitis, sinusitis, stomatitis etc.

